



# i-gel<sup>®</sup> O<sub>2</sub> Resus Pack

Four things in one pack, one less thing to think about



**Airway Management** • Airway Devices



In emergency medicine, you need equipment that's easy, rapid and reliable to use. The i-gel® O<sub>2</sub> Resus Pack contains everything you need to prepare, insert and secure the i-gel® O<sub>2</sub> quickly and efficiently: an i-gel® O<sub>2</sub> supraglottic airway, a sachet of lubricant, airway support strap and a suction tube.


The i-gel® O<sub>2</sub> incorporates a supplementary oxygen port and an integral hook ring for securing of the airway support strap.

The i-gel® O<sub>2</sub> has been designed to facilitate ventilation as part of standard resuscitation protocols, such as those designated by the European Resuscitation Council (ERC) and the American Heart Association (AHA). However, the i-gel® O<sub>2</sub> incorporates a supplementary oxygen port, so can also be used for the delivery of passive oxygenation, or Passive Airway Management (PAM™), as part of an appropriate CardioCerebral Resuscitation (CCR) protocol.



## Airway innovation since 2007

i-gel® is the innovative second generation supraglottic airway device from Intersurgical launched in 2007. The first major development since the laryngeal mask airway, the i-gel® changed the face of airway management and is now widely used in anaesthesia and resuscitation around the world.



### The non-inflatable cuff

The i-gel® gets its name from the innovative soft, gel-like material from which it is made. It is the innovative application of this material that has enabled the development of a unique non-inflatable cuff. This means there is

no need for cuff deflation prior to insertion and no cuff inflation after placement to secure a seal, shortening and simplifying the preparation and insertion procedure.

### Insertion

The i-gel® O<sub>2</sub> is incredibly easy to use. Insertion is rapid and can normally be achieved in less than 5 seconds. The simplicity of i-gel® O<sub>2</sub> means the insertion technique is instinctive, reliable and does not require or benefit from the insertion of fingers or thumbs into the patient's mouth.

### Fixation

The i-gel® O<sub>2</sub> Resus Pack includes a specially designed airway support strap for securing the i-gel® O<sub>2</sub> in position. Fixation is easily achieved by sliding the strap under the patient's neck and securing it in place to the hook ring on the collar of the i-gel® O<sub>2</sub>. This makes it ideal for use where adhesive tape is unsuitable, such as on patients with facial hair.





## i-gel® O<sub>2</sub> Resus Pack

### 15mm connector

Reliable connection to any standard catheter mount or connection

### Proximal end of gastric channel

### Colour coded hook ring

To secure the i-gel® O<sub>2</sub> in position with the airway support strap

### Clearly displayed product information

For quick easy reference. Includes confirmation of size and weight guidance

### Position guide

### Supplementary oxygen port

For the administration of passive oxygenation as a component of cardiocerebral resuscitation (CCR)

### Gastric channel

Enhances patient safety by providing a mechanism for the management of regurgitant fluid

### Integral bite block

Reduces the possibility of airway channel occlusion

### Buccal cavity stabiliser

Aids insertion and eliminates the potential for rotation

### Epiglottic rest

Reduces the possibility of epiglottis 'down folding' and airway obstruction

### Non-inflatable cuff

Eliminates the need for cuff inflation after insertion, allowing easy and rapid insertion

### Distal end of gastric channel



The i-gel® O<sub>2</sub> Resus Pack is available in 3 adult sizes and includes:

- i-gel® O<sub>2</sub> supraglottic airway
- Sachet of lubricant – for quick and easy lubrication of the i-gel® O<sub>2</sub> prior to insertion
- Airway support strap – to secure the i-gel® O<sub>2</sub> in position
- 12 FG suction tube – for insertion through the gastric channel to empty fluid from the stomach.



## Ordering information

The i-gel® O<sub>2</sub> Resus Pack contains an i-gel® O<sub>2</sub>, a sachet of lubricant, airway support strap and suction tube, all packaged in a specially designed clear, sterile, rigid pack. This unique design ensures protection of the contents, both in transit and in storage. The i-gel® O<sub>2</sub> has a colour coded hook ring to allow quick and easy identification of size. The i-gel® O<sub>2</sub> Resus Pack is available in three adult sizes.



Code	Description	Weight	
● AB1015	i-gel® O <sub>2</sub> Resus Pack, large adult – includes a size 5 i-gel® O <sub>2</sub> with orange hook ring, sachet of lubricant, airway support strap and a 12FG suction tube	90+kg	Ⓢ
● AB1014	i-gel® O <sub>2</sub> Resus Pack, medium adult – includes a size 4 i-gel® O <sub>2</sub> with green hook ring, sachet of lubricant, airway support strap and a 12FG suction tube	50–90kg	Ⓢ
● AB1013	i-gel® O <sub>2</sub> Resus Pack, small adult – includes a size 3 i-gel® O <sub>2</sub> with yellow hook ring, sachet of lubricant, airway support strap and a 12FG suction tube	30–60kg	Ⓢ

Ⓢ Sterile

Please note: Suction tube contains PVC.