NIO[™] Adult Quick, Safe, Easy Vascular Access

- Vascular access in under 10 seconds
- Most cost-effective IO solution
- Redundant safety features
- Semi-automatic device
- No exposed needle
- Single-use device
- 5-year shelf life
- Self-contained



Vascular Access in a Flash

When you need immediate vascular access, you can be confident that the NIO[™] Adult is complete and ready to use. The NIO Adult is the quickest, most reliable intraosseous device for gaining immediate vascular access in adult patients.





NIO[™] Adult

Product Highlights

Semi-automatic, self-contained unit facilitates rapid vascular access in under 10 seconds.

Unexposed needle decreases the risk of accidental needle injury.

Redundant safety features prevent accidental activation of the device.

No need for additional accessories makes the NIO Adult a cost-effective solution.

The NIO Adult features an innovative stabilizer hub which serves three functions:

- **Needle stabilizer** assists in securing the needle to the limb
- **Mechanical stopper** limits needle depth, preventing over penetration
- **Stabilizer base** aids in the removal of the catheter

Specifications

Dimensions	139mm (5.47in) L 64.5mm (2.6in) W 43.7mm (1.72in) H
Weight	94.5g
Packaging	Packed individually in a peel-back blister pack
Needle Gauge	15G
Needle Effective Length	42mm (1.65in)
Needle Penetration Depth	25mm (0.98in)
Safety Mechanism	 Rotating safety latch prevents accidental activation Simultaneously compressing the device against the skin while squeezing the trigger wings activates device
Shelf Life*	5 years

Specifications are subject to change without notice. *From date of manufacture



PerSys Medical 5310 Elm St. Houston, TX 77081 USA 1.888.737.7978 info@ps-med.com

www.persysmedical.com