



Australian Government

Department of Health
Therapeutic Goods Administration

Australian Register of Therapeutic Goods Certificate

Issued to

Aero Healthcare

for approval to supply

Aero Healthcare - Thermometer, electronic, intermittent

ARTG Identifier	296501
ARTG Start Date	20/11/2017
Product Category	Medical Device Included Class IIa
GMDN	14035
GMDN Term	Thermometer, electronic, intermittent
Intended Purpose	A hand-held, battery-powered, electronic instrument designed to measure a patient's body temperature. It may comprise an electronic unit with an attached probe or be a single unit (shaped like an ordinary hand-held capillary thermometer) that detects and converts the changes in temperature into variations of some electrical characteristic, e.g., resistance or voltage. These variations of the electrical characteristics are processed in the electronic circuits and in turn displayed, for a short period, as temperature readings. Thereafter the display will automatically turn off or go into standby mode. This is a reusable device

Manufacturer Details	Address	Certificate number(s)
Hangzhou Hua'An Medical & Health Instruments Co Ltd	Building 2 1# Fuzhu Nan RD Wuchang Town Yuhang District , Hangzhou Zhejiang , 310023 China	DV-2017-MC-17604-1

ARTG Standard Conditions

The above Medical Device Included Class IIa has been entered on the Register subject to the following conditions:

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods Act 1989 and Part 5, Division 5.2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Products Covered by This Entry

1. Thermometer, electronic, intermittent

Product Specific Conditions

- No specific conditions have been recorded against this entry.