

QuikClot
**COMBAT
GAUZE**

IT'S **QUIKCLOT** OR IT'S NOT



Control life-threatening bleeding fast
with QuikClot Combat Gauze[®]

TRAUMATIC BLEEDING INJURIES FOR MASSIVE HEMORRHAGE WOUND



FASTER BLEEDING CONTROL

QuikClot® hemostatic products allows you to control bleeding at the scene until medical support arrives

EASY TO USE

Intuitive to use, QuikClot® uses the same placement technique as standard gauze

UPDATED PACKAGING

Compact, cube-sized packaging, and eight tear points expedite access to dressing for quicker application

ACTIVATES

Kaolin Activates Factor XII (Intrinsic Pathway)

ACCELERATES

Accelerates the Clotting Cascade

FORMATION

Robust Clot Formation

CONTROL

Faster Bleeding Control

Aero Healthcare Pty Ltd

63 Seaton Street, Armidale, NSW, 2350

T 61 2 6776 7200

F 61 2 6771 3369

E sales@aerohealthcare.com

WWW.AEROHEALTHCARE.COM

