



TECHNICAL FEATURES

Product size (L x W x H): 183 x 45 x 5cm

Net weight: 7.5kg

Gross weight: 9kg

Load bearing capacity: Less than 159kg

DESCRIPTION

Made of durable PE material with no discharge contaminator, this versatile spine board is equipped with safe belts and can double as a floatable device for water rescues. Its completely translucent structure allows for X-rays, and it's compatible with most head immobilization devices and strap mechanisms.

APPLIED FOR

This versatile spine board is ideal for use in emergency situations in hospitals, sports facilities, ambulances, and outdoors. Its durable and flexible construction, along with compatibility with most head immobilization devices and strap mechanisms, makes it easy to transport and secure patients. It can also be used as a floatable device for water rescues.

CAUTIONS:

1. Do not exceed the load bearing capacity of less than 159kg.
2. Do not use the spine board if it is damaged or shows signs of wear and tear.

MAINTENANCE

1. Keep clean as routine (incl. sterilization).
2. Often check if any parts loose or not.

STORAGE & TRANSPORTATION

1. Store this product in the place of damp proof & non-corrosion.
2. Common transporting vehicle can carry this product.