

Australian Government

Department of Health and Aged Care

Therapeutic Goods Administration

Public Summary

Summary for ARTG Entry: 156690 Stryker Australia Pty Ltd - Defibrillator, semi-automated

ARTG entry for Medical Device Included Class IIb

Sponsor Stryker Australia Pty Ltd

Postal Address PO Box 970, ARTARMON, NSW, 1570

Australia

ARTG Start Date 12/11/2008

Product Category Medical Device Class IIb

Status Active

Approval Area Medical Devices

Conditions

- The inclusion of the kind of device in the ARTG is subject to compliance with all conditions placed or imposed on the ARTG entry. Refer Part 4-5, Division 2 (Conditions) of the Therapeutic Goods (Medical Devices) Regulations 2002 for relevant information.
- Breaching conditions of the inclusion related to the device of the kind may lead to suspension or cancellation of the ARTG entry; may be a criminal offence; and civil penalties may apply.

Manufacturers

Name Address

HeartSine Technologies Ltd 207 Airport Road

West Belfast, Northern Ireland, BT3 9ED

United Kingdom

Products

1. Defibrillator, semi-automated

Product Type Single Device Product Effective Date 12/11/2008

GMDN 37805 Defibrillator, semi-automated

Intended Purpose A device that can analyse an electrocardiogram (ECG) to determine whether or not a defibrillator shock should be delivered. It is

attached to the patient through a pair of adhesive defibrillation electrodes that serve both to monitor the ECG and to deliver the

defibrillator discharges. The semi-automated defibrillator indicates to the operator when to deliver the shock.

Specific Conditions

No Specific Conditions included on Record

© Commonwealth of Australia. This work is copyright. You are not permitted to re-transmit, distribute or commercialise the material without obtaining prior written approval from the Commonwealth. Further details can be found at http://www.tga.gov.au/about/website-copyright.htm.