TRANSPARENT PLASTERS




AERO PLAS[™]

TRANSPARENT

The AeroPlast^M Transparent range brings to you a transparent aesthetic film, blending with every skin tone and providing your wound with a waterproof barrier in a plaster form with a full range of sizes.





PROTECTIVE WRAPPER

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening.

Individually wrapped for guaranteed quality and tamper protection.

FEATURES


- Transparent
- Waterproof
- Conforming
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond[™] Adhesive Technology





TRANSPARENT PLASTERS PRODUCT INDICATIONS

The AeroPlast[™] Transparent range is available with a selection of shapes and sizes to ensure the perfect fit in different applications.



A TRANSPARENT AESTHETIC FILM PLASTER BLENDS WITH ANY SKIN TONES FOR A SUBTLE APPLICATION.

AEROPLAST[™] TRANSPARENT LAYER

The AeroPlast[™] transparent film layer offers an aesthetic wound covering that subtly blends with all skin colours. This layer provides conforming, waterproof protection against dirt and germs.

SUPABOND® ADHESIVE

Exclusively engineered for unrivaled latex free and hypoallergenic adhesion.

WOUND PAD

Super absorbent non-woven wound pad that helps reduce skin maceration & aid in faster healing.

NON STICK BARRIER

A non-adherent barrier, provides wound protection and minimises trauma on removal.

PROTECTIVE FILM

PE coated & silicone release papers ensure easy removal without touching the wound pad & helps maintain superior adhesion.

HEALTHCARE | WE MAKE IT BETTER