




VISUAL PLASTERS

The AeroPlast™ Visual Plaster range is a premium latex free, hypo-allergenic, visually detectable woundcare range specifically designed for the retail food industry and non-commercial food preparation.



FEATURES

- Visually detectable
- Food grade
- Waterproof
- Conforming
- Breathable
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond[™] adhesive technology
- Ideal for the retail food industry



PROTECTIVE WRAPPER

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening.

Individually wrapped for guaranteed quality and tamper protection.







VISUALLY DETECTABLE PLASTERS PRODUCT INDICATIONS

The AeroPlast™ Visual plasters standard size is larger for even greater wound isolation. The Visual range also offers a dressing strip that can be cut to size to suit larger or awkward wounds.

X-WIDE STRIP 2.5 X 7.5 cm

Actual size
Represented at 100% when printed to A4

Actual size


Actual size

Actual size

Actual size

Actual size

THE ONLY COLOUR THAT DOES NOT OCCUR NATURALLY IN FOOD IS BLUE!



AEROPLAST™ VISUAL LAYER

A vibrant blue plastic layer that bends, stretches and conforms to suit all applications. Providing waterproof protection against dirt and germs while maintaining breathability for the wound.

SUPABOND® ADHESIVE

Exclusively engineered for unrivaled latex free and hypoallergenic adhesion.

WOUND PAD

Super absorbent non-woven wound pad that helps reduce skin maceration & aid in faster healing.

NON STICK BARRIER

A non-adherent barrier, provides wound protection and minimises trauma on removal.

PROTECTIVE FILM

PE coated & silicone release papers ensure easy removal without touching the wound pad & helps maintain superior adhesion.