


**AERO** PLAST™

**PLASTIC**

AeroPlast™ Plastic offer an economical, washproof range of conforming bandages perfect for protecting wounds from dirt and germs that can cause infection.



**PROTECTIVE WRAPPER**

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening.

Individually wrapped for guaranteed quality and tamper protection.

**FEATURES**


- Washproof
- Conforming
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond™ Adhesive Technology




# PLASTIC PLASTERS

## PRODUCT INDICATIONS

The AeroPlast™ Plastic range is available with a selection of shapes and sizes to ensure the perfect fit in different applications.



**THE ECONOMICAL, WASHPROOF, SHEER AESTHETIC DRESSING  
RECOMMENDED FOR ALL GENERAL APPLICATIONS.**



### AEROPLAST™ PLASTIC LAYER

A unique plastic layer that bends, stretches and conforms to different shapes and applications. A washproof layer that protects the wound from dirt and germs, while being air permeable allowing it to breath.

### SUPABOND® ADHESIVE

Exclusively engineered for unrivaled latex free and hypoallergenic adhesion.

### WOUND PAD

Super absorbent non-woven wound pad that helps reduce skin maceration & aid in faster healing.

### NON STICK BARRIER

A non-adherent barrier, provides wound protection and minimises trauma on removal.

### PROTECTIVE FILM

PE coated & silicone release papers ensure easy removal without touching the wound pad & helps maintain superior adhesion.