


**AERO** PLAST™


**PREMIUM METAL DETECTABLE**

AeroPlast™ Premium Detectable bandages are easily identifiable with their vibrant blue colour. Manufactured with a unique blue material chosen for its vibrant colour and exceptional waterproof qualities and conformability. Designed to withstand resistance to water, oil and grease and with its metal detectable qualities, this bandage is perfect for food processing industries.



**FEATURES**

- Metal detectable
- Visually detectable
- Waterproof
- Ultra conforming
- Breathable
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond™ adhesive technology
- Ideal for the manufacturing & processing sector



**PROTECTIVE WRAPPER**



A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening.

Individually wrapped for guaranteed quality and tamper protection.




# PREMIUM METAL DETECTABLE PLASTERS PRODUCT INDICATIONS

The AeroPlast™ Premium Detectable range is available with a selection of shapes and sizes to ensure the perfect fit in different applications.



**A LONGER LASTING, MORE DURABLE PLASTER THAT CAN ACTUALLY DELIVER MEASURABLE REDUCTIONS IN PLASTER USAGE.**



### AEROPLAST™ VISUALLY DETECTABLE LAYER

A layer of premium material, chosen for its unique qualities that outperform others in essential areas of stretch, conformability and comfort. Its ability to resist against water, oil and grease, ensures less chance of lifting, allowing it to perform for longer. The layer completely seals the wound from bacteria, while maintaining breathability.

### SUPABOND® ADHESIVE

Exclusively engineered for unrivaled latex free and hypoallergenic adhesion.

### METAL DETECTABLE FOIL

A metal detectable wound pad backing enables guaranteed detection where metal detection equipment is operational, making this plaster ideal for the manufacturing & processing sector and the ultimate high security.

### WOUND PAD

Super absorbent non-woven wound pad that helps reduce skin maceration & aid in faster healing.

### NON STICK BARRIER

A non-adherent barrier, provides wound protection and minimises trauma on removal.

### PROTECTIVE FILM

PE coated & silicone release papers ensure easy removal without touching the wound pad & helps maintain superior adhesion.