

POWERED BY



# BLEEDING CONTROL CATALOGUE

Aero Healthcare's comprehensive bleeding control pricelist  
including tourniquets, haemostatic dressings & bleed kits



**A TRAUMA  
SOLUTION ISN'T  
SUPPLEMENTARY  
TO FIRST AID,  
IT'S ESSENTIAL.**

## CONTENTS

*	4	TOURNIQUET SPECIFICATIONS
*	5	TOURNIQUETS
*	6-7	BLEEDING CONTROL KIT SMALL
*	8-9	BLEEDING CONTROL KIT MEDIUM
*	10	CONSUMABLES
*	11	ISRAELI BANDAGES
*	12	TRAINING
*	13	BAGS



INTUITIVE. STRONG. LIFESAVING.

# DESIGNED TO WORK, **FAST.**

The RapidStop® Tourniquet is an award-winning, rapid application tourniquet that maximises the chances of survival. Developed at MIT and perfected by EMS, military and industry specialists. RapidStop® represents a haemostatic solution with military-grade ruggedness, all age suitability, superior single-handed application and the fastest average occlusion time.



# SPECIFICATION SHEET

Ratchet mechanical advantage of 8:1; Easier and quicker to achieve high compression forces, particularly necessary for application on upper leg (thigh)

Simple, rapid application around trapped limbs: with just gross motor skills, the RapidStop® Tourniquet Quick-Clasp can be detached and reattached for application around a trapped limb - no dexterous strap routing required. 'QUICK-CLASP' can be opened and closed easily even when wearing gloves

Intuitive product-integrated application instruction on Label (Step 1)



TIME RECORD LABEL LOCATED BEHIND D-RING STRAP, BELOW QUICK-CLASP

Instantaneous one-handed cinch vastly improves time to occlusion in all application scenarios. Easy and quick to apply; requires no stabilizing hand and supplies a **mechanical advantage of 2:1** prior to effecting tightening mechanism [NOTE: a poorly-performed initial cinch is a typical cause for a tourniquet application failure.]

Large, rugged D-ring 5.1cm / 2" width (3.8cm / 1.5" height) makes it simple and quick to grasp



**LEVER** Intuitive product-integrated application instruction on Lever (Step 2)

**RATCHET**

**STRAP**

**QUICK-CLASP**

**LABEL**

**D-RING**

Minimum 600N tensile strength: 200N is typically the largest tensile force needed in a tourniquet to fully occlude blood flow in large thigh application; RapidStop® Tourniquet is designed with a minimum 3x factor of safety

RapidStop® Tourniquet is designed for limb sizes circumference from 17.8cm / 7" (83.8cm / 33") this equates to the 1st percentile of female upper arm sizes up to the 99th percentile of male upper thigh sizes; encompassing the vast majority of adult limb sizes

Mil-W-17337 spec; 3.8cm / 1.5" wide with 8,000N breaking strength and superior abrasion resistance

TYPE	MILITARY	FIRST RESPONDER	TRAINING
PART#	RST100	RST101	RST102
COLOR	BLACK	ORANGE	BLUE
UPC#	9341394017504	9341394017511	9341394017528
TEMPERATURE RANGE	-50c to + 50c (-58f to +122f)		
TENSILE STRENGTH	minimum 600N		
PACKAGE SIZE	19.56cm long x 5.1cm wide x 3.8cm high (7.7" long x 2" wide x 1.5" high)		
WEIGHT	127.6g (4.5oz)		
STRAP	Mil-W-17337 ; 3.81cm width ; 8,000N		
LIMB CIRCUMFERENCE	17.8cm to 83.8cm (7" to 33")		
RESISTANCE	IMPACT-RESISTANT TO LARGE CRUSHING, BENDING FORCES		

# TOURNIQUETS

RAPIDSTOP®

01

## INTUITIVE.

Fast and easy single-handed self-application are paramount to a tourniquet's effectiveness in any environment. RapidStop® stands out from the competition as one of the few tourniquets that exceed these life-saving requirements utilising key patented and proprietary technologies.

02

## STRONG.

Made in America and built tough to military standards. The patented ratchet lever technology is made from the same polycarbonate used in the F-22 jet fighter canopy. The tourniquet's band is made from high-strength nylon webbing rated to a two-ton breaking point. RapidStop® has constructed the tourniquet to endure extreme temperatures and resist any corrosion, ensuring it stands true in the toughest of conditions.

03

## LIFESAVING.

Easy as 1-2-3. Three quick and simple steps set RapidStop® Tourniquet apart from the competition. Clear instructions enable an untrained user to apply the tourniquet with speed and reliability. Integrated instructions on the RapidStop® Tourniquet are stitched on as weather-resistant waterproof vinyl tags. Pocket-sized instruction cards make learning and becoming familiar with the RapidStop® easy.



## BLACK TOURNIQUET

MILITARY | RST100 | \$52.50, 25pcs x \$45.15  
RRP \$72.00  
Carton QTY: 25 | Min Order: 1



## ORANGE TOURNIQUET

FIRST-RESPONDER | RST101 | \$52.50, 25pcs x \$45.15  
RRP \$72.00  
Carton QTY: 25 | Min Order: 1



## BLUE TOURNIQUET

TRAINING | RST102 | \$52.50, 25pcs x \$45.15  
RRP \$72.00  
Carton QTY: 25 | Min Order: 1

# BLEEDING CONTROL KIT

SMALL

BASED ON THE GLOBALLY ACCEPTED CONTENTS FOR A PERSONAL CARRY, THE SMALL RAPIDSTOP BLEEDING CONTROL KIT IS DESIGNED TO GIVE THE USER CAPABILITIES TO ATTEND TO A PATIENT WITH MAJOR HAEMORRHAGING. INCLUDED IN THE CONTENTS IS THE WORLD-RENOWNED ISRAELI BANDAGE AND THE LATEST IN TOURNIQUET TECHNOLOGY FROM RAPIDSTOP.

## CONTENTS

	QTY
RAPIDSTOP® Tourniquet Orange Strap	1
FIRSTCARE™ Civilian Trauma & Hemorrhage Control Bandage 10 cm x 17 cm	1
AEROINSTRUMENT™ Stainless Steel Universal Shears 19 cm	1
AEROGLOVE™ Large Nitrile Powder-Free Gloves Pair	1
AEROSUPPLIES™ Marker Pen	1



Designed for quick and effective control of mass haemorrhaging, RapidStop Bleeding Control Kit provides unmatched occlusion times thanks to the RapidStop Tourniquet.

The 5.11 Tactical softpack bags are ideal for use as a stand-alone kit, and the retail-ready packs can be utilised as a bleed module or a refill kit. These kits are available in small and medium sizes.



### QUICK ACCESS PACK

RSK401 | \$63.00 | RRP \$90.00

MIN ORDER: 4

The RapidStop Bleeding Control Pack is designed for a rapid-access bleed kit with clear identifying markings and clear backing that allows contents to be viewed. This pack can be used as a stand-alone kit, a bleed module or a refill pack.



### TACTICAL BAG

RSK410 | \$141.75 | RRP \$190.00

MIN ORDER: 2

The medium tactical bleed kit is packed out into the UCR IFAK pouch from 5.11 Tactical®. Single-hand access design & drop down front pockets make the pouch ideal for a rapid access bleed kit. Internal pockets have bungee cord tie-downs to keep equipment secure and drain holes make the kit easily washable if it becomes dirty or contaminated with fluids.

The MOLLE attachment system is utilised on the front and back for easy customisation.

# BLEEDING CONTROL KIT

MEDIUM

DESIGNED TO PROVIDE THE USER WITH THE VITAL TOOLS TO SAVE A LIFE, THE MEDIUM KIT CONTAINS THE ITEMS REQUIRED TO ATTEND TO MAJOR HAEMORRHAGE AND SUCKING CHEST WOUNDS. THIS KIT GIVES THE RESPONDER THE MEANS TO STABILISE THE PATIENT FOR THE 'GOLDEN HOUR'.

## CONTENTS

	QTY
RAPISTOP® Tourniquet Orange Strap	1
FIRSTCARE™ Civilian Trauma & Hemorrhage Control Bandage 10 cm x 17 cm	1
HYFIN® Vent Chest Seal – Twin Pack	1
AEROINSTRUMENT™ Stainless Steel Universal Shears 19 cm	1
AERORESCUE™ Thermal Blanket 140 cm x 210 cm	1
AEROGLOVE™ Large Nitrile Powder-Free Gloves Pair	2
AEROSHIELD™ Disposable Face Shield	1
AEROSUPPLIES™ Marker Pen	1



Designed for quick and effective control of mass haemorrhaging, RapidStop Bleeding Control Kit provides unmatched occlusion times thanks to the RapidStop Tourniquet.

The 5.11 Tactical softpack bags are ideal for use as a stand-alone kit, and the retail-ready packs can be utilised as a bleed module or a refill kit. These kits are available in small and medium sizes.



### QUICK ACCESS PACK

RSK411 | \$99.75 | RRP \$140.00

MIN ORDER: 4

The RapidStop Bleeding Control Pack is designed for a rapid-access bleed kit with clear identifying markings and clear backing that allows contents to be viewed. This pack can be used as a stand-alone kit, a bleed module or a refill pack.



### TACTICAL BAG

RSK410 | \$141.75 | RRP \$190.00

MIN ORDER: 2

The medium tactical bleed kit is packed out into the UCR IFAK pouch from 5.11 Tactical®. Single-hand access design & drop down front pockets make the pouch ideal for a rapid access bleed kit. Internal pockets have bungee cord tie-downs to keep equipment secure and drain holes make the kit easily washable if it becomes dirty or contaminated with fluids.

The MOLLE attachment system is utilised on the front and back for easy customisation.

# CONSUMABLES



## QUIKCLOT® COMBAT GAUZE

QC-CG | \$83.00 | RRP: \$115.50

Min Order: 1



## QUIKCLOT® EMS 4X4 DRESSING

QC-4X4D | \$24.94 | RRP: \$40.00

Min Order: 1



## QUIKCLOT® EMS ROLL

QC-EMSR | \$44.50 | RRP: \$60.00

Min Order: 1



## HYFIN® VENT CHEST SEAL

TWIN PACK | ARHYV1 | \$36.47 | RRP: \$45.00

Min Order: 1



## BLACK TACTICAL GLOVES

PAIR | AGNBK02 | \$1.26 | RRP: \$2.00

Min Order: 50



## PACKING GAUZE

RSH340 | \$5.68 | RRP: \$7.35

Min Order: 6



## BLIZZARD CASUALTY BLANKET

BH-01 | \$231.94 | RRP: \$399.99

Min Order: 1



## BLIZZARD EMS BLANKET

BPS-210 | \$77.01 | RRP: \$149.95

Min Order: 1



## BLIZZARD EMS BABY WRAP

BPS-23 | \$22.58 | RRP: \$44.95

Min Order: 1

# ISRAELI BANDAGES



## 4" MILITARY WOUND DRESSING

FCP-01 | \$12.88 | RRP: \$22.99

10 x 17cm | Min Order: 12



## 4" CIVILIAN WOUND DRESSING

FCP-03 | \$12.88 | RRP: \$22.99

10 x 17cm | Min Order: 12



## 6" MILITARY WOUND DRESSING

FCP-02 | \$16.56 | RRP: \$29.55

15 x 18cm | Min Order: 12



## 6" CIVILIAN WOUND DRESSING

FCP-06 | \$16.56 | RRP: \$29.49

15 x 18cm | Min Order: 12



## 8" MILITARY WOUND DRESSING

MULTI TRAUMA | FCP-09T | \$32.66 | RRP: \$57.89

30 x 30cm | Min Order: 12



## 8" CIVILIAN WOUND DRESSING

ABDOMINAL | FCP-10+ | \$27.16 | RRP: \$47.99

30 x 30cm | Min Order: 12



## 4" TACTICAL WOUND DRESSING

MILITARY TRAUMA TREATMENT | FCP-T3 | \$16.60 | RRP: \$29.45

10 x 17cm | Min Order: 12

# TRAINING



## RAPIDSTOP® TOURNIQUET

TRAINING | RST102 | \$52.50 | RRP: \$72.00

Carton QTY: 25 | Min Order: 1



## PACKING GAUZE

TRAINING | RSH340 | \$5.68 | RRP: \$7.35

Min Order: 6



## HAEMORRHAGE ARM

DUAL CUT | HCS01 | \$2,058.00 | RRP: \$2,225.00

Min Order: 1



## TOURNIQUET TRAINER

SINGLE ARM | ETT001A | \$1,538.25 | RRP: \$1,660.00

Min Order: 1



## TOURNIQUET TRAINER

SINGLE LEG | ETT001L | \$1,538.25 | RRP: \$1,660.00

Min Order: 1



## STABILISING BRACKET

FITS HCS01 | HCS01A | \$157.50 | RRP: \$170.00

Min Order: 1



## LACERATION TRAINER

LIMB | HCT03 | \$756.00 | RRP: \$815.00

Min Order: 1



## LACERATION TRAINER

BASIC SQUARE | HCT04 | \$525.00 | RRP: \$565.00

Min Order: 1



## GUNSHOT WOUND TRAINER

BASIC SQUARE | HCT04A | \$525.00 | RRP: \$565.00

Min Order: 1

# BAGS



## 5.11 UCR IFAK POUCH

MEDIUM - BLACK | 511-UIP | \$49.35 | RRP \$49.95

Min Order: 1



## TOURNIQUET POUCH

FITS RAPIDSTOP | RST105 | \$18.90 | RRP \$28.00

Min Order: 1



## 5.11 FLEX MED POUCH

BLACK | 511-FMP | \$52.49 | RRP \$68.75

Min Order: 1



## 5.11 FLEX TACMED POUCH

BLACK | 511-FTP | \$55.13 | RRP \$69.90

Min Order: 1



## 5.11 UCR THIGH RIG

BLACK | 511-UTR | \$57.75 | RRP \$59.99

Min Order: 1



## 5.11 UCR SLING PACK

BLACK | 511-USP | \$147.00 | RRP \$179.99

Min Order: 1

POWERED BY



**AERO HEALTHCARE AU PTY LTD**

63 SEATON STREET, ARMIDALE NSW 2350 AUSTRALIA  
PO BOX 1635, ARMIDALE NSW 2350 AUSTRALIA  
FREECALL 1800 628 881  
TEL +61 2 6776 7200 | FAX +61 2 6771 3369  
SALES@AEROHEALTHCARE.COM | WWW.AEROHEALTHCARE.COM

HEALTHCARE | WE MAKE IT BETTER

FIRST AID CONSUMABLES

PRE-EMERGENCY

FIRST AID KITS

WOUND CARE