

# Hybrid Trauma Kit

Aero Healthcare's exclusive Hybrid Trauma Kit is designed to combine Defibrillator & First Aid tools to equip rescuers with essential items to treat major bleeds to the arms and legs including sudden cardiac arrest following a traumatic event.

Hard Case Hybrid Trauma Kit for HeartSine Defibrillators ensures your defibrillator is protected against extreme weather, whilst remote working.



Burns Module

Trauma Module

Eye Wound Module

Dressing and Bandage Module

Your choice of HeartSine Samaritan Pad with 8-year Warranty

**Hard Case** ■ Tested to IP67 ■ MIL-STD-810F for Vibrations and Impact Resistance ■ Crushproof ■ Dust-proof ■ Waterproof ■ Lightweight ■ Bespoke Thermal Protection Foam

Applying the correct wound care is the first step in effective first aid, the Hybrid Trauma Kit makes it easy!



HeartSine Samaritan PAD with CPR Advisor model 500P.

- ✓ **Built-in confidence.** We recommend the SamaritanPAD 500P uses an Impedance Cardiogram (ICG) to assess the effectiveness of the CPR and give the rescuer feedback on compressions.
- ✓ **User-friendly.** Easy-to-understand visual and voice prompts guide a user through the entire process.
- ✓ **Always ready.** A System Status Ready Indicator flashes to show that the complete system is operational and ready for use. The AED also automatically runs a regular self-check.

The comprehensive wound care is stocked in an easy-to-use package with colour and clear labels to ensure the right first aid is applied for the right injuries.



## MODULE KIT CONTENTS

Aero Tuffcut Scissors 7.5x190mm - each  
AeroPore Clear Adhesive Tape 2.5x5cm - each  
AeroBurn Sterile Gel Dressing 20x20cm - each  
AeroBurn Sterile Face Gel Dressing Face Mask - each  
Aero Security Seal - Pack 10  
AeroHazard Clinical Waste Bag 381x711x990mm - each  
Rapidstop Tourniquet - Orange Strap  
Aero A6 First Aid Leaflet - each  
AeroWash Eyewash 250ml - each  
AeroPlast w/p Assorted - Env10  
AeroGlove Nitrile Gloves Powder free extra large - pair  
AeroShield CPR Face Shield - each  
AeroPlast Emergency Foil Blanket silver 140x210cm - each  
AeroWound Trauma Dressing Medium 10x18cm pad - each  
AeroBand CC Cotton Calico Tri-Bandage 90x90x127 - each  
AeroWound Trauma Dressing Large 15x18cm pad - each  
AeroGlove Nitrile Gloves Powder (free medium) - pair  
AeroGlove Nitrile Gloves Powder (free large pair) - pair  
AeroWipe Saline Cleansing Wipes (bulk) - each



HEALTHCARE | WE MAKE IT BETTER

CALL **01403 790 704**

EMAIL **sales@aerohealthcare.co.uk**

VISIT **aerohealthcare.com/uk**

© Aero Healthcare