



# AGNPF01E PR (PAIR)

Blue nitrile powder free disposable glove



### Sectors:

- Automotive • Construction & Utility • Food Manufacturing and Preparation • Janitorial • Oil and Gas • Medical and Dental

### Features

**Manufactured from:** Nitrile butadiene

**Latex free:** reducing the risk of hypersensitivity and allergic reactions

**Low in residual chemicals:** To reduce the risk of contact dermatitis associated with these chemicals

**Flexible material:** For improved feel and fit reducing fatigue during extended wear

**Strength:** Offers better puncture, tear and abrasion resistance than natural rubber or vinyl

**Powder free:** Reduces the risk of dust contamination

**Textured finish:** For improved grip

**Food Contact:** Tested in accordance with European standard EN1186 making it suitable for food contact

**Colour**  
**Sizes**

**Blue**  
**S-XL**



Size	Product Code
S/7	AGNPF01E-S-PR
M/8	AGNPF01E-M-PR
L/9	AGNPF01E-L-PR
XL/10	AGNPF01E-XL-PR

### Standards

Contains accelerators which may cause allergic reactions

Manufactured in accordance with Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425

Compliant with the Medical Devices Regulation (EU) 2017/745 Class 1

Tested to European Standards EN455: Medical Gloves for Single Use parts 1,2,3&4

Product Food Approved to EC1935/2004

EU Type Examination carried out by SATRA

Tested to European Standards EN420:2003+A1:2009 & EN ISO 374-1:2016 & EN ISO 374-5:2016

### Packing

1 pair of gloves per pack



Disposable Gloves

Nitrile Disposables