

BOLIN CHEST SEAL

- ✓ MDR Compliant
- ✓ UKCA Marked & CE Marked



The Bolin Chest Seal is a sterile occlusive chest wound dressing for treating open pneumothorax and preventing tension pneumothorax that results from gunshots, stab wounds, or other penetrating chest trauma. The patented (US Patent No. 7,834,231) failure-proof triple-valve design of the BCS allows air and blood to escape while preventing the re-entry of either, thereby eliminating any unwanted gas or liquid exchange at the trauma site.

SPECIFICATION:



Dimensions: 7.5" x 8" x .25"



Weight: 2.1 oz



Product shelf life: 5 years

AVAILABLE IN:



- ✓ Failure-proof triple-valve design
- ✓ Rugged polyurethane disc
- ✓ Covered with a thick layer of jell-based adhesive