




CR2 Replacement Battery

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The Lifepak CR2 battery replacement is fully compatible with all Lifepak CR2 defibrillators, including both semi-automatic and fully automatic models, ensuring that your device remains operational in critical situations. It is recommended to replace the battery every four years or when prompted by the device's maintenance commands, to ensure the defibrillator is always ready for immediate use. The battery is designed for easy and quick replacement, minimizing downtime and ensuring that the defibrillator can continue to deliver its reliable and fast heart rhythm analysis. Lifepak CR2 defibrillators are known for their exceptional performance, offering vital features such as CPR coaching and metronome sounds to assist first responders in providing effective chest compressions. The straightforward and user-friendly process of changing the battery further enhances the defibrillator's readiness, ensuring it remains a dependable tool in emergency situations.

PHE1014



- ✓ Replace every 4 years
- ✓ Fast Heart Rhythm analysis
- ✓ Simple battery replacement process

For this Defibrillator:



Lifepak
CR2