

INFRARED THERMOMETER

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



The AeroDiagnostic™ Infrared Thermometer ensures accurate body temperature readings. This device can measure temperature when positioned 5 centimetres away from someone's forehead. Used for temperature management, it allows operators to measure someone's temperature without the need for physical contact. An accurate and stable reading is delivered within one second, offering a speedy response time. The AeroDiagnostic™ Infrared Thermometer supports medical and healthcare professionals that need to document patient body temperatures, allowing them to note temperatures without invasion of personal space. More accurate and timely readings can allow medical professionals and first aiders to make quick decisions in regard to patient safety and care. Likewise, for those looking to document the temperatures of visitors on worksites in response to COVID-19 or other illnesses, the device ensures simple and easy readings within seconds of use.

AVAILABLE IN:



- ✓ 12 Months
- ✓ Non-Contact
- ✓ 1 second response time
- ✓ 5cm measuring distance