

# OxyMulti-Mask™

High Oxygen • High Humidity • Aerosol Therapy



## Patient Safety

- Reduces the risk of CO<sub>2</sub> rebreathing
- Reduces the risk of emesis aspiration
- For nose or mouth breather

## Patient Comfort

- Reduces feelings of claustrophobia
- No build up of heat or humidity under the mask
- Unrestricted communication
- Access for fluids and oral medication

## Oxygen Therapy Compliance

Uninterrupted oxygen therapy aids compliance of prescribed oxygen therapy

## Operational Efficiencies

One OxyMask™ replaces multiple traditional oxygen therapy interfaces:  
SKU standardization • Inventory reduction • Reduced medical waste  
• Reduces complexity for bedside staff

# OxyMulti-Mask™

- Oxygen flows from 1 to 15+ litres per minute - flush
- Delivers 23% to 83% FiO<sub>2</sub>
- Delivers oxygen through a standard small bore tube, or using a 22mm OD adapter attached to a universal swivel elbow
- The universal elbow can be connected to most aerosol nebulizers to deliver effective aerosol therapy and not compromise O<sub>2</sub> delivery
- Can be used with, or without high humidity or aerosol therapy
- Ideal for acute care/EMS situations when you need high FiO<sub>2</sub>, high humidity and/or delivery of medicated aerosol therapy
- Not made with natural rubber latex or phthalates such as DEHP
- 25 masks/case



OxyMulti-Mask™  
Adult



Code	Description	Mask Size	Tubing Length	Elastic
OHH-1025-8	OxyMulti-Mask™ with adapter & 7' tubing	Adult	7'	Single