

AED RESCUE RUCKSACK

- ✓ UKCA Marked
- ✓ CE Marked



The AED Rescue Rucksack is designed to carry any HeartSine® Samaritan® public access defibrillator (PAD) and its accessories to any location. This robust rucksack enables the convenient transport of a life-saving heart-safe solution, ensuring quick AED access during sudden cardiac arrest emergencies. Features the international AED symbol and the word 'AED' for easy identification.

- 600D/420D water resistant polyester.
- Wide zip opening defibrillator compartment.
- Padded adjustable shoulder straps with attachment clip.
- Easy grab zip pullers for ease of opening in an emergency situation.
- Breathable padded back panel.
- Elasticated bottle side pocket.
- Metal ventilation eyelets.
- Padded grab handle.
- External dimensions: (H) 50 cm x (W) 30 cm x (D) 18 cm. Internal capacity: 23 litres.

