

NEEDLE PICK UP KIT 1APP REFILL, 0.3L SHARPS BIN

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



This compact sharps bin has been designed with sustainability, environmental impact, ease of use, efficiency and compliance in mind. Manufactured from 40% recycled plastic and certified to ISO23907:2019 for sharps injury protection, UN3291 medical or clinical wastes as well as compliant with HTM 07-01 guidance for the safe management of healthcare waste.

The flip top cover ensures easy and safe disposal featuring a non-return petal opening ensuring any disposed needles or other sharp items are safely contained. It accepts up to 5ml syringes. Comes pre-assembled with essential medical disposal accessories included. Ideal for healthcare, lab, dental, and other industries that require safe disposal of sharps. Sharps bin measures 159mm (H) x 79mm (L) x 40mm (W).

Kit Contents:

- | | |
|--|---------------------------|
| 1 x 0.3 Litre Sharps Disposal Bin | 1 x Pair Plastic Tweezers |
| 1 x Yellow Self-Seal Waste Bag (355mm x 380mm) | 1 x Instruction Card |
| 1 x Bioguard Hand & Surface Wipe | |



SK-1APPR

