

## 4-WHEEL PLUS PATIENT FOLDING TRANSIT WHEELCHAIR

### FOR UNEVEN GROUND

- ✓ MDD Compliant
- ✓ UKCA Marked & CE Marked



4-Wheel Plus Transit Wheelchair is an enhanced version of the basic Transit Wheelchair for transferring a person requiring additional comfort and support across rough and uneven surfaces or carry patients up or down a flight of stairs. A minimum of two operators are required to carry the chair but additional may be necessary depending on the weights to be lifted. Because of the lifting and Manual Handling involved, it is recommended that full Risk Assessment and Personal Emergency Evacuation Plan (PEEP) should be carried out prior to purchase.

#### SPECIFICATION:

- Chair Weight 13kg
- Chair Length 90cm
- Chair Width 58cm
- Chair Height 128cm
- Load capacity 130kg



#### AVAILABLE IN:



- ✓ Large rear wheels enable movement across rough and uneven surfaces.
- ✓ Acts like a wheelchair for easy manoeuvrability when on level ground.
- ✓ Flat seat, arms, and footrest help to relax the client
- ✓ Integrated handles ensure the chair is carried in the correct position