

TRUECLOT®

IPT Inguinal Packing Trainer

The TrueClot Inguinal Packing Trainer is a wearable training aid that simulates a junctional hemorrhage that cannot be treated with a tourniquet. The trainer can be placed on a manikin or worn by an actor for a higher level of realism to teach advanced wound packing skills.

Features

Wearable, life-threatening wound simulating an inguinal hemorrhage untreatable with a standard tourniquet

Fully packable with a 12 ft strip of simulated hemostatic dressing

Compatible with junctional tourniquets (when used on a manikin)

Easy to clean with water and air dry

Durable ballistic nylon construction

Includes 1L of TrueClot Blood Simulant, 5 packs of gauze, and storage case

Available in Dark and Light Flesh tones



MADE IN THE USA

Dimensions

Kit in case: 11 x 13.5 x 5.2 in (28 x 34.3 x 13.2 cm) - 3.9 lbs (1.8 kg)

Wound trainer: 10.5 x 5.25 in (15.9 x 22.2 cm) - Wound depth: 2.5 in (6.4 cm)

Liquid capacity: 33.8 oz (1 liter)

Item # 5174

www.trueclot.com