


AERO PLAST™

TRANSPARENT

The AeroPlast™ Transparent range brings to you a transparent aesthetic film, blending with every skin tone and providing your wound with a waterproof barrier in a plaster form with a full range of sizes.



PROTECTIVE WRAPPER

A cold-sealed paper wrapper ensures a sterile bandage, with distinct opening tabs for easy fumble-free opening.

Individually wrapped for guaranteed quality and tamper protection.

FEATURES

- Transparent
- Waterproof
- Conforming
- Sterile
- Latex Free
- Hypoallergenic
- SupaBond™ Adhesive Technology




TRANSPARENT PLASTERS

PRODUCT INDICATIONS

The AeroPlast™ Transparent range is available with a selection of shapes and sizes to ensure the perfect fit in different applications.



**A TRANSPARENT AESTHETIC FILM PLASTER
BLENDS WITH ANY SKIN TONES FOR A SUBTLE APPLICATION.**



AEROPLAST™ TRANSPARENT LAYER

The AeroPlast™ transparent film layer offers an aesthetic wound covering that subtly blends with all skin colours. This layer provides conforming, waterproof protection against dirt and germs.

SUPABOND® ADHESIVE

Exclusively engineered for unrivaled latex free and hypoallergenic adhesion.

WOUND PAD

Super absorbent non-woven wound pad that helps reduce skin maceration & aid in faster healing.

NON STICK BARRIER

A non-adherent barrier, provides wound protection and minimises trauma on removal.

PROTECTIVE FILM

PE coated & silicone release papers ensure easy removal without touching the wound pad & helps maintain superior adhesion.