



Data Sheet

ABOUT FIRST CARE PRODUCTS

First Care Products provides innovative first-response medical products which save lives by improving field hemorrhage control, and also yield significant savings in time and resources. The company's unique, combat and clinically proven direct pressure Emergency Bandage has been adopted by military and civilian organizations worldwide.

ULTIMATE FIELD HEMORRHAGE CONTROL

- Combat and clinically proven direct pressure bandage
- Trauma bandage and secondary dressing all in one
- Quick, easy and effective application; and self-application, even one-handed

SAVING LIVES WITH A UNIQUE DIRECT PRESSURE SYSTEM

According to military statistics, more than half of all military combat deaths result from blood loss. Now, The Emergency Bandage raises combat hemorrhage control to the next level, consolidating numerous pieces of treatment equipment into a single effective and easily applied unit.

- 'Standard' in major military organizations worldwide

Learn more at: www.aerohealthcare.com/uk

FEATURES

- Integral pressure applicator exerts immediate direct pressure to the wound
- Secondary sterile dressing maintains pad and pressure firmly in place and lowers risk of further infection
- Built-in closure bar; no pins, clips, tape, velcro, or knots
- Non-adherent pad eliminates pain and wound re-opening upon removal
- Compact; durable package

BENEFITS

- All-in-one bandage lowers the amount of equipment military medical personnel need to carry
- Easily applied by untrained personnel - ideal personal field dressing for every soldier
- Unique Life Cycle Cost (LCC) administrative, logistical and training savings compared to existing multiple-bandage treatment solutions - no need to order, stock, supply, and train on numerous pieces of equipment
- Rapid control of bleeding stabilizes wounded and saves lives

PRODUCT SPECIFICATIONS



64950

Description. Military personal field dressing

Pad size. 4 inch / 10 cm

Color. Green




64960

Description. Military personal field dressing

Pad size. 6 inch / 15 cm

Color. Green



Learn more at: www.aerohealthcare.com/uk